

## Introduction

Reversed-phase HPLC is an invaluable tool also for the analytical and preparative separation of peptides and proteins. Owing to the availability of different pore sizes and particle sizes, the alkyl-bonded silica gel products are economically the first choice for both analytical and preparative separations.

In purification of peptides and proteins, dynamic axial compression system is the most useful method to purify the compound. Using the DAC column, it is easily fractionated by changing the separation conditions. (e.g. packing material, column length, column volume, etc) We reported that choosing optimum pore size of gel gave good peak shape and separation in analytical scale. Based on the results, we applied the analytical conditions to preparative/process scale separation. Effective separation was obtained similarly to the analytical separation of peptides and proteins.

1

## Optimum Packing Conditions for 50 mml.D. DAU Column

### Packing conditions

Packing material : YMC \* GEL ODS-A (15 $\mu$ m, 120 Å)  
Weight of packing material : 250 g  
Column size : 200 x 50 mml.D.  
Packing pressure : 6.4 MPa  
Slurry solvent : 100% methanol  
Concentration of slurry : 35%


### HPLC conditions for inspection of packing procedure

Eluent : methanol / water (85/15)  
Flow rate : 50 mL/min  
Temperature : ambient  
Detection : UV at 254 nm  
Sample : 1. Toluene  
2. Methyl benzoate

5

## Impact of Pore Size on Efficiency (1)

### Angiotensin II (MW 1,046) in analytical




120 Å pore size is most efficient at all the loading levels.

9

## Loadability (1)

### Insulin (MW 5,700) in preparative

YMC \* GEL C8 (15 $\mu$ m, 200 Å)



Up to 400 mg loaded amount, peak shapes are sharp as the analytical separation.

13

## Preparative Separation of Insulin (2)

### Loadability of insulin

Loaded amount / mg	Purity / %	Recovered amount / mg	Recovery / %
500	99.2	441	88.9
700	99.3	561	80.7

Choosing the optimum gel, high purity and high yield of recovered insulin was obtained.  
Effective separation was obtained similarly to the analytical separation. It could be scaled up to gram scale purification.


17

## New Dynamic Axial Compression Column - DAU series -

- Easily packing and unpacking via inlet and outlet tube
- Unit is available for high (70 to 100 bar) pressure
- Column diameter extends to 600 mm
- Compact design with built in packing station
- Dynamic axial bed compression yields densely packed beds
- Recommended media: RPC media

2

## Reproducibility of Repacking the Gel




By using a 50 mml.D. DAU column, repacking procedure was attempted. After 7 times repacking, plates/m and Tf are still good as the initial state.

6

## Impact of Pore Size on Efficiency (2)

### Insulin (MW 5,700) in analytical




200 Å pore size is the best choice for samples up to the 0.1 mg loading level.

10

## Loadability (2)

### Ovalbumin (MW 45,000) in preparative

YMC \* GEL C4 (15 $\mu$ m, 300 Å)



At all loading levels, high efficiency was obtained.


14

## Conclusions

- After 7 times repacking procedure, peak shapes and column performance are as good as the initial state.
- In preparative scale, it is also important to choose the optimum pore size and the right ligand to achieve optimal separation of peptides or proteins, similarly to analytical separation.
- You can use wide variation of gels from YMC Co. Ltd. [e.g. three pore sizes (120, 200, 300 Å), several particle sizes, several ligand groups] We also offer some dimension of dynamic axial compression column.

18


## DAU Products Line



3

## Mechanical Stability of the Gel

### Column pressure




After 7 times repacking, column pressure is almost same as initial state. It would be no formation of fines. This shows the gel is stable under the flow rate.  
High mechanical stability ensures longer lifetime of the gel.

7

## Impact of Pore Size on Efficiency (3)

### Ovalbumin (MW 45,000) in analytical




300 Å pore size is most efficient at all the loading levels.  
At low loading levels, 200 Å pore size also shows good efficiency.

11

## Scale Up from Analytical to Preparative


YMC \* GEL C8, 200 Å



Each samples are separated in preparative conditions, similarly to analytical scale.

15

## Packing Steps for a DAU Column



Easily packing and unpacking via inlet and outlet tube

4

## Peptide and Protein Purification by Reversed Phase Silica Gel

Small organic molecules are retained/eluted by a distribution mechanism. On the other hand, peptides and proteins are retained/eluted by an adsorption-desorption (on-off) mechanism. Due to this mechanism, the pore size plays a key role in determination of resolution and loading amount in separation of peptides and proteins.

Based on the results in analytical separation as we reported previously, we attempted to separate peptides and proteins in preparative conditions. This study shows scalability of the separation using DAU column.

8

## Optimized Stationary Phase for Separation


MW		C18	C8	C4	
5,000	120 Å	⊙	○	△	⊙ : excellent ○ : good △ : moderate
20,000	200 Å	○	⊙	○	
100,000	300 Å	△	○	⊙	

C18 column with 120 Å pore size is suitable for small peptides up to MW 5,000, similar to the analyses of ordinary small molecules. In the case of large peptides or small proteins up to MW 20,000, the C8 column with 200 Å pore size gives the best efficiency. Most of proteins are eluted effectively by C4 column with 300 Å.

12

## Preparative Separation of Insulin (1)

YMC \* GEL C8, 200 Å



16

## New Dynamic Axial Compression Column

DAU-450 with slurry container

