

Jeffrey A, Kakaley^{1*}, Masatoshi Taniguchi², Chiyuki Awahara², Takuya Sagami², Takashi Sato² and Noritaka Kuroda²
¹YMC America, Inc, ²YMC CO., LTD.

Introduction

Antibody drugs are expected to enhance the effect of therapies without adverse side effects. During the manufacturing processes of immunoglobulin G (IgG), which is a main ingredient of antibody drugs, its aggregates and fragments could be generated as impurities. Size exclusion chromatography (SEC) analysis is known as a common method for determination of IgG characteristics as well as quality control in the production processes of IgG. YMC has developed a new SEC column with 2 µm particle size and Diol bonded phase. High resolution analysis of aggregates and fragments derived was performed with this new SEC column. Furthermore, smooth method transfer from conventional 5 µm SEC column to this new 2 µm SEC column could also be achieved while reducing analysis time by using a shorter column. In this poster, we will introduce these advantages of this 2 µm Diol SEC column through several examples of IgG aggregates and fragments analysis.


Basic characteristics

Specifications of novel Size Exclusion Diol column

| Product name | Pore size (Å) | Particle size (µm) | Base material | Functional group | Pressure limit |
|--------------|---------------|--------------------|---------------|------------------|---------------------------|
| Diol-200 | 200 | 2 | Porous silica | Dihydroxy propyl | 2 µm: 6530 psi (45 MPa) |
| | | 3 | | | 3,5 µm: 2900 psi (20 MPa) |
| Diol-300 | 300 | 5 | | | |

- Silica gel based size exclusion chromatography (SEC) columns bonded with dihydroxypropyl group.
- High pressure tolerance that is sufficient to withstand operating pressure required for high throughput SEC analyses
- Very low nonspecific adsorption of biomacromolecules, including monoclonal antibodies

Calibration curves and separation characteristics



| # | Sample | MW | # | Sample | MW |
|---|---------------|---------|----|--------------------|--------|
| 1 | IgM | 900,000 | 9 | α1-Antrypsin | 50,000 |
| 2 | Thyroglobulin | 670,000 | 10 | Ovalbumin | 45,000 |
| 3 | IgA | 390,000 | 11 | Carbonic anhydrase | 30,000 |
| 4 | Fibrinogen | 340,000 | 12 | Trypsin inhibitor | 20,100 |
| 5 | γ-Globulin | 158,000 | 13 | Myoglobin | 17,000 |
| 6 | IgG | 150,000 | 14 | α-Lactalbumin | 14,100 |
| 7 | Transferrin | 75,000 | 15 | Ribonuclease A | 13,700 |
| 8 | HSA | 66,000 | 16 | Cytochrome c | 12,400 |


Column : YMC-Pack Diol-200 / Diol-300, 300 X 4.6 mm i. d.
 Eluent : 0.1 M KH₂PO₄-K₂HPO₄ (pH 7.0) containing 0.2 M NaCl
 Flow rate : 0.165 ml/min
 Temp. : 25°C
 Detection : UV at 280 nm

- If pore size is the same, calibration curves were identical across particle sizes (2, 3 and 5 µm). This feature would realize seamless method transfer between particles (HPLC=UHPLC).
- As expected, 300 Å pore material gave higher resolution on high molecular weight side, while 200 Å pore material had advantage on separation of low molecular weight side. 300 Å pore material is suitable for antibody aggregates analysis, and 200 Å pore material works best for antibody fragment analysis.


Analysis of monoclonal antibody and its aggregates

Influence of flow rate on antibody analysis


Diol-300, 2 µm, 300 X 4.6 mm i. d.



Improvement of resolution / analysis throughput by using the 2 µm column



Evaluation of column durability




- Stability of 2 µm YMC-Pack Diol column on monoclonal antibody analysis was evaluated. Performance of the column was maintained for 500 injections, and expected to last much longer.

- Theoretical plate count and resolution of monoclonal antibody increased as flow rate was decreased. This suggests that lower flow rate is preferable if scope is high resolution separation analysis. On the other hand, higher flow rate is suitable for increasing analysis throughput, especially when resolution is sufficient.
- The backpressure at 0.8 ml/min with 300 mm length by 4.6 mm i. d. column was within the pressure limit.

- 2 µm YMC-Pack Diol in 300 X 4.6 mm i. d. column greatly improved resolution between aggregates and monomer peak compared to 5 µm YMC-Pack Diol in the same dimension.
- On 2 µm, 150 mm length column, resolution between antibody monomer and aggregates was the same as the 5 µm, 300 mm length column. This would suggest that analysis time can be reduced by half just by changing the column from 5 µm, 300 mm length column to 2 µm, 150 mm length column.
- Operating pressure was low enough to run with conventional HPLC system. *Special care is needed in terms of minimizing system volume in order to avoid extra-column band spreading when using at low flow rate.

Lot-to-Lot reproducibility




- Separation reproducibility of 2 µm YMC-Pack Diol column between manufacturing lots was evaluated. Separation was identical between lots. This feature would contribute to robust method development.

Analysis of monoclonal antibody fragments

Method of Papain digestion

- 50 mM KH₂PO₄-K₂HPO₄ (pH 7.5) buffer was prepared.
- Papain and other components were added to the buffer at the following concentrations.
 - Papain 0.1 g/l
 - NaCl 0.15 mol/l
 - L-cystein 0.01 mol/l
 - EDTA 0.01 mol/l
- A 3 mg/mL sample of monoclonal antibody was made using solution 2 as the diluent.
- The solution was stored for 17.5 hours at 25 °C.

Overlaid chromatogram of the intact and the digested monoclonal antibody



- The analysis of papain digested monoclonal antibody was performed with the novel 2 µm YMC-Pack Diol-200 column. Intact antibody and fragments were well separated.

Conclusions

The new 2 µm Diol SEC columns exhibited improved resolution and allowed higher throughput of monoclonal antibody analyses.

This feature would offer significant advantages for antibody analyses focused on aggregate and fragment control in research as well as quality control of pharmaceutical products.