

Introduction

Nucleic acids such as antisense, siRNA and aptamer are expected as next-generation pharmaceuticals following antibody drugs. For providing these drugs, purification and separation analysis that can recognize slight structural differences after synthesis are important issues.

Non-porous anion exchange column is generally suitable for analysis of oligonucleotides. Thus we tried to optimize an analysis method of single-stranded DNA and RNA of about 20 mer, using BioPro IEX QF which is a nonporous high performance anion exchange column. For optimization, we changed some conditions such as type of mobile phase and column temperature. As the result, good separations could be obtained for the oligonucleotides with single-base difference in length.


In this poster, we will introduce further details about the optimization of separation conditions.

Specifications of BioPro IEX columns

	BioPro IEX QF	BioPro IEX SF	BioPro IEX QA	BioPro IEX SP
Matrix	non-porous hydrophilic polymer beads		porous hydrophilic polymer beads	
Particle size (µm)	3, 5		5	
Charged group	-CH ₂ N ⁺ (CH ₃) ₃	-CH ₂ CH ₂ CH ₂ SO ₃ ⁻	-CH ₂ N ⁺ (CH ₃) ₃	-CH ₂ CH ₂ CH ₂ SO ₃ ⁻
Counter ion	Cl ⁻	Na ⁺	Cl ⁻	Na ⁺
Ion-exchange capacity (meq/ml-resin)	0.075-0.110	0.230-0.290	0.075-0.100	0.070-0.095
Dynamic binding capacity (mg/ml-resin)	>12 (BSA)	>10 (human-IgG)	>110 (BSA)	>70 (human-IgG)
Usable pH range	2-12			
Column size (length X i.d.(mm))	30 X 4.6, 50 X 4.6, 100 X 4.6			

Sample Group 1 (Phosphodiester oligonucleotides ; PO)


1	Single-stranded DNA (ssDNA)	5'-TCATCACAAGTACCAAT-3' (DNA 20 mer)
2	Single-stranded DNA (ssDNA)	5'-GTCATCACAAGTACCAAT-3' (DNA 21 mer)
3	Single-stranded RNA (ssRNA)	5'-UCAUCACACUGAAUACCAU-3' (RNA 20 mer)
4	Single-stranded RNA (ssRNA)	5'-UCAUCACACUGAAUACCAU-3' (RNA 21 mer)
5	Single-stranded RNA (ssRNA)	5'-U(M)C(M)A(M)U(M)C(M)A(M)C(M)A(M)C(M)U(M)G(M)A(M)A(M)U(M)A(M)C(M)C(M)A(M)A(M)U(M)-3' (2'-Ome RNA 20 mer)
6	Single-stranded RNA (ssRNA)	5'-G(M)U(M)C(M)A(M)U(M)C(M)A(M)C(M)A(M)C(M)U(M)G(M)A(M)A(M)U(M)A(M)C(M)C(M)A(M)A(M)U(M)-3' (2'-Ome RNA 21 mer)



N(M)=2'-Ome RNA

Sample Group 2 (Phosphorothioate oligonucleotides ; PS)

10	Single-stranded DNA	5'-TATATATATATATATATATATATTT-3' (DNA 15 mer 12PS, 2PO)
11	Single-stranded DNA	5'-TATATATATATATATATATATATTT-3' (DNA 15 mer 13PS, 1PO)
12	Single-stranded DNA	5'-TATATATATATATATATATATATAT-3' (DNA 15 mer All PS)
13	Single-stranded RNA	5'-U [^] C [^] A [^] A [^] U [^] C [^] A [^] C [^] A [^] C [^] A [^] U [^] G [^] A [^] A [^] U [^] A [^] C [^] A [^] A [^] U [^] -3' (RNA 20 mer All PS)
14	Single-stranded RNA	5'-G [^] A [^] C [^] A [^] U [^] C [^] A [^] C [^] A [^] C [^] U [^] G [^] A [^] A [^] U [^] A [^] C [^] A [^] A [^] U [^] -3' (RNA 21 mer All PS)




^=Phosphorothioated

(i) Improvement of carryover peak

If initial gradient concentration of NaCl was low(ex. 50 mM), carryover occurred. But increased initial gradient concentration of NaCl up to 400 mM enabled to avoid carryover with good reproducibility.


(ii) Improvement of peak tailing




(i ~ ii) Summary

By using BioPro IEX QF, the condition for the analysis of oligonucleotides (PO) was optimized. We concluded 10 mM NaOH(as buffer solution), NaClO₄(as salt) and higher initial gradient concentration of salt were preferable to suppress carryover and peak tailing.


(iii) Examples of analysis under the optimized condition




(iv) Separation of ssDNAs with single-base difference in length (differ in type of base of 5'end of DNA 21mer)




(v) Resolution of phosphorothioate oligonucleotides with different degrees of thiolation



(vi) Difference in required salt concentrations for eluting modified RNA (All PS) and normal RNA (All PO)



(vii) Separation optimization trial of phosphorothioate oligonucleotides with single-base difference in length



Conclusions

- By using BioPro IEX QF:
 - Each ssDNA, ssRNA and 2'-Ome ssRNA with single-base difference in length can be successfully separated.
 - [All PS], [13PS, 1PO] and [12PS, 2PO], which consist of 15-mer ssDNA, can be also separated under the optimized condition.
 - Higher salt concentration is required to elute All PS compared to eluting All PO.
 - All PS with single-base difference in length can be separated to some extent. So we continue to establish more optimized conditions in future.