

## 1. Introduction

Recently, ultra-fast LC has been widely applied to analysis of various fields such as pharmaceuticals, foods, and environmental. Although commercial columns packed with sub-2  $\mu\text{m}$  or 2  $\mu\text{m}$  particles give improved efficiency and reduced analysis time, acceptable retention or separation has not been often achieved in many hydrophilic compounds such as metabolites or natural compounds.

To solve these problems, we have lined up a new category of products, YMC-UltraHT Hydrosphere C18, into our 2  $\mu\text{m}$  columns for ultra-fast LC. The balanced hydrophilic/hydrophobic nature of this packing material provides strong retention and superior selectivity of highly polar compounds. Also this new column maintains all the advantages of 2  $\mu\text{m}$  YMC-UltraHT Pro C18 column. Our 2  $\mu\text{m}$  columns show almost same efficiency of sub-2  $\mu\text{m}$  columns with about 40 % lower back-pressure and they can be used with ordinary LC systems.

In this poster, we will show some example cases of fast and efficient separation of pharmaceuticals and metabolites containing polar compounds. We will also compare the retention and elution achieved with Hydrosphere C18 and those achieved with other C18 columns for ultra-fast LC.

## 2. Features of YMC-UltraHT columns

### High column efficiency with minimum column pressure

- Functionalized ultra-pure 2.0  $\mu\text{m}$  silica gel with 120 Å pore size.
- Designed for lower pressure but better performance than sub-2  $\mu\text{m}$  particles.
- Instruments are not only specific ultra-fast LC but ordinary LC available.

### Easy scalability of separation

- Identical selectivity as 3  $\mu\text{m}$  or 5  $\mu\text{m}$  conventional Pro series columns within same bonded phase.
- Applicable for various compounds such as pharmaceuticals, foods and natural products.
- Simple method transfer applicable without changing eluent conditions.

## 3. New UltraHT column for polar compounds

### YMC-UltraHT Hydrosphere C18


- Designed for strong retention and superior resolution selectivity of polar compounds.
- Excellent reproducibility of retention time under 100 % aqueous condition.
- Unique selectivity that differs from standard ODS, Pro C18.

	Hydrosphere C18	Pro C18
feature	for polar compounds	standard type ODS
particle size ( $\mu\text{m}$ )	2, 3, 5	2, 3, 5, 10
carbon contents (%)	12 %	16 %
bonded phase	monomeric	monomeric
end-capping	completely end-capped	completely end-capped


## 4. Fast and efficient separation of polar compounds on YMC-UltraHT Hydrosphere C18

### Ultra-fast analysis under 100 % aqueous condition

#### Xanthine oxidase (XOD) inhibitor and related metabolites



#### Easy method transfer from conventional LC to Ultra-fast LC




#### Comparison of retention stability in 100% aqueous mobile phase

##### Evaluation method

Analytical condition  
same as the optimized method in above figure

##### Test procedure

- analysis at 0.6 ml/min (initial: drawn in blue)
- decrease to 0.2 ml/min and flow continuously for 30 min (with lower pressure)
- analysis at a 0.6 ml/min (after decreasing flow: drawn in red)
- stop the pump for 30 min (with no pressure)
- analysis at 0.6 ml/min (after stopping pump: drawn in green)




- In both cases of decreasing flow rate and stopping the pump, the considerable loss of retention is observed with Brand G2 (standard ODS) and Brand I11 (ODS recommended to use in 100% aqueous mobile phase). YMC-UltraHT Hydrosphere C18 shows excellent stability and reproducibility.

- It is important to choose appropriate column for the reproducible analysis in a highly aqueous mobile phase, because 0.2 ml/min is common flow rate for 2.0 or 2.1 mm i. d. column.


#### Superior selectivity comparing to commercial sub-2 $\mu\text{m}$ columns for polar compounds

sample: spiked human serum




### Ultra-high resolution of polar metabolites

#### Metabolites of organic solvents in urine




#### Improved resolution of positional isomers



- By changing particle size from 5  $\mu\text{m}$  to 2  $\mu\text{m}$ , the resolution between each peak is significantly improved.
- Furthermore, the 100 mm length column with higher efficiency can achieve baseline separation of *p*- and *m*-isomer of methylhippuric acid (peak 6 and peak 7) which are hardly separable with conventional reversed-phase columns.

#### Analysis of human urine sample with 2 $\mu\text{m}$ 100 X 2.0 mm i. d. column



## 5. Conclusions

- YMC-UltraHT Hydrosphere C18 shows same selectivity as existing 3  $\mu\text{m}$  and 5  $\mu\text{m}$  Hydrosphere C18, so easy method transfer can be achieved between conventional LC and Ultra-fast LC without changing eluent condition.
- YMC-UltraHT Hydrosphere C18 provides superior resolution and excellent retention reproducibility of highly polar compounds even in 100 % aqueous mobile phase.
- The combination of high efficiency and unique selectivity of YMC-UltraHT Hydrosphere C18 enables the high-throughput analyses of various pharmaceuticals, foods and natural products containing polar compounds.