



YMC-Pack Pro C18




pH-stable 2–8
ideal for basic, acidic and polar compounds

YMC-Pack Pro C18 RS



pH-stable 1–10
hydrophobic compounds
steric selectivity

Hydrosphere C18




pH-stable 2–8
also 100% aqueous eluents


Easy ProFamily Column Selection

for Method Development

YMC Radar Chart




YMC-Pack Pro C8



pH-stable 2.0–7.5
good alternative to C18 phases

YMC-Pack Pro C4



pH-stable 2.0–7.5
ideal for faster separations



For more information about the YMC-Pack Pro series, scan the QR code or visit:

<https://www.ymc-america.com/ymc-brand/ymc-pack-pro/>



8 Charlestown Street, Devens, MA 01434
Tel: 978.487.1100 • info@ymc-america.com
www.ymc-america.com